Introduction

What would happen to your cells if they made a poisonous chemical? You might think that they would die. In fact, your cells are always making poisonous chemicals. They do not die because your cells use enzymes to break down these poisonous chemicals into harmless substances. Enzymes are proteins that speed the rate of reactions that would otherwise happen more slowly. The enzyme is not altered in the reaction. You have hundreds of different enzymes in your cells. Each of these is responsible for one particular reaction that occurs in the cell.

In this lab, you will study an enzyme that is found in many living tissues. The name of the enzyme is **catalase** (KAT-uhlayss); it speeds up a reaction which breaks down hydrogen peroxide, a toxic chemical, into 2 harmless substances: water & oxygen. The reaction is as follows:

$$2H_2O_2 \rightarrow 2H_2O + O_2$$

This reaction is important to cells because **hydrogen peroxide** (H_2O_2) is produced as a byproduct of many normal cellular reactions. If the cells did not break down the hydrogen peroxide, they would be poisoned & die.

In this lab, you will study the catalase found in liver cells. You will be using chicken or beef liver that your teacher purchased in the supermarket. It might seem strange to use dead cells to study the function of enzymes. This is possible because when a cell dies, the enzyme remains intact & active for several weeks, as long as the tissue is kept refrigerated. Recall that the substrate is the molecule that the enzyme acts on, and the products are the molecules produced by the reaction.

Answer the pre-lab questions before moving on to Part A:

Pre-lab Ouestions:

- 1. What is the name of the enzyme studied in this lab?_
- 2. What chemical does the catalase break down? (In other words: what is the substrate?)
- 3. Catalase breaks the hydrogen peroxide into what two harmless substances? (In other words: what are the products?)

Catalase Lab Question #1: Are enzymes reusable?


Procedure:

• Place 20 mL (about the height of your thumbnail) of 3% hydrogen peroxide solution into a clean test tube & add a small piece of liver.

1. Observe the bubbles. What gas is being release?

Throughout this investigation you will estimate the rate of the reaction (how rapidly the solution bubbles) on a scale of 0-5 (0 = no reaction, 1 = slow, ...5 = very fast). Assume that the reaction that you just saw above proceeded at a rate of "4".

Recall that a reaction that absorbs heat is endothermic (reactants have less energy; products have more which will feel hot); A reaction that gives off heat (reactants have more energy; products have less & feel cold) is exothermic.



Type of reaction: energy-absorbing reaction

Type of reaction: energy-releasing

Now, feel the temperature of the test tube with your hand.

- 2. Has the tube gotten warmer or colder? _____ Is the reaction endothermic or exothermic? _____
 - Pour off the liquid in the tube into the sink.

the state of the s	4. Was there another reacting hydrogen peroxide to the reward why?	emaining liver?			
Catalase is present isn many kinds of living tissue. You will now test for the presence of catalase in tissues other than liver. Place 20 ml of hydrogen peroxide in each of 3 clean test tubes and then add each of the three test substances to the tubes. As you add each test substance, record the reaction rate (0-5) for each tube. Substance Rate of Reaction (0-5) Potato Apple or Carrot Chicken 1. Which tissues contained catalase? 2. Do some contain more catalase than others? How can you tell? Catalase Lab Question #3: How are enzymes in living tissue affected by temperature? Under certain conditions enzymes are denatured. An enzyme is denatured when the protein molecule loses its proper shape & cannot function. Some things that can denature an enzyme are high temperatures, extremes of pH, heavy metals & alcohol. → Put a piece of liver into the bottom of a clean test tube & cover it with a small amount of distilled water. Place this test tube in a boiling water bath (100°C) for 5 minutes. What will boiling do to an enzyme? • Remove the test tube from the hot water bath, allow it to air cool, then pour out the water. Add 20 mL of hydrogen peroxide. CAUTION: use a test-tube holder when handling hot test tubes. 2. What is happening in the test tube? • Record the reaction rate in the first column of the data chart (0 = no reaction → 5 = high reaction). Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) O'C (freezing) Rate of Reaction (0→5) • Put equal quantities of liver into 3 clean test tubes and a 10 mL. H ₂ O ₂ into 3 other test tubes. Put one test tube of liver & one of H ₂ O ₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest). 4. Why did the reaction proceed slowly at 0°C?	5. Are enzymes reusable?				
liver. Place 20 ml of hydrogen peroxide in each of 3 clean test tubes and then add each of the three test substances to the tubes. As you add each test substance, record the reaction rate (0-5) for each tube. Substance Potato Apple or Carrot Chicken 1. Which tissues contained catalase? 2. Do some contain more catalase than others?	Catalase Lab Questions #	#2: What Tissues Con	ntain Catalase?		
Apple or Carrot Chicken 1. Which tissues contained catalase? 2. Do some contain more catalase than others?	liver. Place 20 ml of hydro	gen peroxide in each	of 3 clean test tubes and the	hen add each of the three	
Apple or Carrot Chicken 1. Which tissues contained catalase? 2. Do some contain more catalase than others?	Substance		Rate	of Reaction (0-5)	
1. Which tissues contained catalase? 2. Do some contain more catalase than others? How can you tell?					
1. Which tissues contained catalase? 2. Do some contain more catalase than others?		t			
2. Do some contain more catalase than others?	Chicken				
Catalase Lab Question #3: How are enzymes in living tissue affected by temperature? Under certain conditions enzymes are denatured. An enzyme is denatured when the protein molecule loses its proper shape & cannot function. Some things that can denature an enzyme are high temperatures, extremes of pH, heavy metals & alcohol. → Put a piece of liver into the bottom of a clean test tube & cover it with a small amount of distilled water. Place this test tube in a boiling water bath (100°C) for 5 minutes. 1. What will boiling do to an enzyme? • Remove the test tube from the hot water bath, allow it to air cool, then pour out the water. Add 20 mL of hydrogen peroxide. CAUTION: use a test-tube holder when handling hot test tubes. 2. What is happening in the test tube? • Record the reaction rate in the first column of the data chart (0 = no reaction → 5 = high reaction). Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) • Put equal quantities of liver into 3 clean test tubes and a 10 mL H ₂ O ₂ into 3 other test tubes. Put one test tube of liver & one of H ₂ O ₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H ₂ O ₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest). 4. Why did the reaction proceed slowly at 0°C? 4. Why did the reaction proceed slowly at 0°C?	1. Which tissues contained	catalase?			
Under certain conditions enzymes are denatured. An enzyme is denatured when the protein molecule loses its proper shape & cannot function. Some things that can denature an enzyme are high temperatures, extremes of pH, heavy metals & alcohol. → Put a piece of liver into the bottom of a clean test tube & cover it with a small amount of distilled water. Place this test tube in a boiling water bath (100°C) for 5 minutes. 1. What will boiling do to an enzyme? • Remove the test tube from the hot water bath, allow it to air cool, then pour out the water. Add 20 mL of hydrogen peroxide. CAUTION: use a test-tube holder when handling hot test tubes. 2. What is happening in the test tube? • Record the reaction rate in the first column of the data chart (0 = no reaction → 5 = high reaction). Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) 0°C (freezing) e Put equal quantities of liver into 3 clean test tubes and a 10 mL H₂O₂ into 3 other test tubes. Put one test tube of liver & one of H₂O₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H₂O₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest). 4. Why did the reaction proceed slowly at 0°C?	2. Do some contain more of	catalase than others? _	How can you tell?		
Under certain conditions enzymes are denatured. An enzyme is denatured when the protein molecule loses its proper shape & cannot function. Some things that can denature an enzyme are high temperatures, extremes of pH, heavy metals & alcohol. → Put a piece of liver into the bottom of a clean test tube & cover it with a small amount of distilled water. Place this test tube in a boiling water bath (100°C) for 5 minutes. 1. What will boiling do to an enzyme? • Remove the test tube from the hot water bath, allow it to air cool, then pour out the water. Add 20 mL of hydrogen peroxide. CAUTION: use a test-tube holder when handling hot test tubes. 2. What is happening in the test tube? • Record the reaction rate in the first column of the data chart (0 = no reaction → 5 = high reaction). Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) 0°C (freezing) e Put equal quantities of liver into 3 clean test tubes and a 10 mL H₂O₂ into 3 other test tubes. Put one test tube of liver & one of H₂O₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H₂O₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest). 4. Why did the reaction proceed slowly at 0°C?	Catalase Lab Ouestion #	3: How are enzymes	in living tissue affected b	ov temperature?	
 Record the reaction rate in the first column of the data chart (0 = no reaction → 5 = high reaction). Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature	 → Put a piece of liver test tube in a boiling 1. What will boiling do to • Remove the test tu hydrogen peroxide 	ng water bath (100°C) an enzyme? be from the hot water c. CAUTION: use a te	bath, allow it to air cool, st-tube holder when hand	then pour out the water.	
Data Table: Rates of Enzyme Reaction at Different Temperatures: Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) 0°C (freezing) Rate of Reaction (0→5) • Put equal quantities of liver into 3 clean test tubes and a 10 mL H ₂ O ₂ into 3 other test tubes. Put one test tube of liver & one of H ₂ O ₂ into each of the following water baths: • Ice bath (0°C) • Room temperature bath (about 22°C) • Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H ₂ O ₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest). 4. Why did the reaction proceed slowly at 0°C? 4. Why did the reaction proceed slowly at 0°C? 4. Why did the reaction proceed slowly at 0°C?				a montion \rightarrow 5 — high m	anation)
Temperature 100°C (boiling) 37°C (body temp) 22°C (room temp) 0°C (freezing) • Put equal quantities of liver into 3 clean test tubes and a 10 mL H ₂ O ₂ into 3 other test tubes. Put one test tube of liver & one of H ₂ O ₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H ₂ O ₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest)				$\frac{10 \text{ reaction } \rightarrow 3 = \text{ mgn re}}{2}$	eaction).
 Put equal quantities of liver into 3 clean test tubes and a 10 mL H₂O₂ into 3 other test tubes. Put one test tube of liver & one of H₂O₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H₂O₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest)	Temperature			22°C (room temp)	0°C (freezing)
tubes. Put one test tube of liver & one of H ₂ O ₂ into each of the following water baths: ✓ Ice bath (0°C) ✓ Room temperature bath (about 22°C) ✓ Warm water bath (37°C) (human body temperature) • After 3 minutes, pour each tube of H ₂ O ₂ into the corresponding tube of liver & observe the reaction. Record the rates of reaction in the last three columns of the data chart (0→5). 3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest)		o of liver into 2 aloon	test tubes and a 10 ml. II	O into 2 other test	
3. What is the optimum temperature for catalase? (This is the temperature at which the reaction proceeds fastest) 4. Why did the reaction proceed slowly at 0°C?	tubes. Put one test ✓ Ic ✓ Ro ✓ W • After 3 minutes, p observe the reaction	tube of liver & one of e bath (0°C) com temperature bath arm water bath (37°C) our each tube of H_2O_2	(about 22°C) (human body temperaturinto the corresponding tu	llowing water baths: re) lbe of liver &	A A A
	3. What is the optimum ter			at which the	•
	4. Why did the reaction pro	oceed slowly at 0°C?_			
		•			
	5. Why did the reaction no	t proceed at all at 100	· C !		